

TORAY


Hollow Fiber Dialyzer

Toraysulfone **TS-S/U** SERIES



Toraysulfone—Excellent Polysulfone Membrane

Optimum Membrane Structure by Advanced Nano-Technology



Sharp Molecular Weight Cut-off

Polysulfone membrane has an asymmetric membrane structure with high solute removability and ultrafiltration.

Fig. 1 shows, A: TEM observation on membrane section, B and C: AFM observation on membrane surface. All of them show the asymmetric pore structure of the membrane.

Table 1 is a comparison of each polysulfone membrane structure by nanoscopic characterization using Atomic Force Microscope.

The result suggests that Toray Toraysulfone has the following characteristics compared to other polysulfone membranes.

1. Substances can pass easily through the membrane because of its thin skin layer, small tortuosity and short pore length.
2. Larger molecular weight substances such as albumin do not easily pass through Toraysulfone membrane because of its low water content and small pore at skin layer.

Toray polysulfone "Toraysulfone" has an optimal membrane structure among high-flux polysulfone membranes.

Crosslinked Structure of PVP in Toraysulfone Membrane

Polyvinylpyrrolidone (PVP) in Toraysulfone membrane is crosslinked during γ -ray sterilization, and less PVP is eluted from membrane.

Fig. 2 is an observation of three different polysulfone membranes after soaking in DMAc (solvent for polysulfone).


PVP dissolves away together with polysulfone in A2 and A3 membrane. However crosslinked PVP in Toraysulfone remains as a transparent structure.

Crosslinked PVP in Toraysulfone remains even in the solvent.

Toraysulfone has less elution of PVP compared to other membranes.

Table 1

Membrane	Water content (%)		Tortuosity		Thickness (µm)		Pore diameter (nm)		Pore length (µm)
	Skin layer	Support layer	Skin layer	Support layer	Skin layer	Support layer	Skin layer	Support layer	
Toraysulfone	27	70	1.13	1	2	38	8.8	418	40.3
A1	31	70	1.73	1	3	42	11	494	47.2
A2	42	73	1.14	1	7	38	9.5	499	46.0
A3	47	77	1.80	1	2	38	13.0	699	41.6





Spacer Yarns


Toraysulfone hollow fibers are covered by "Spacer yarns" as shown in Fig. 3.

Spacer yarns facilitate the dialysate to flow uniformly around the hollow fibers and to reduce the "mass transfer resistance in dialysate side" without affecting a pressure drop in the dialysate compartment.

Spacer yarns help to enhance the efficiency of dialysis.



New design




Helical Slits Structure

Toraysulfone TS series has specially designed baffle structure with helical slits, in order to obtain the uniform flow of dialysate. (Fig. 4)

In this new design, the baffle was arranged so as to surround both ends of the fiber bundle. The dialysate uniformly penetrates into the fiber bundle from surrounding slits.

Uniformity of dialysate flow in the new Toraysulfone TS series can be observed by X-ray CT scan as shown in Fig. 5.



S
Clin
rent
pati
C
cant
In
stud
clea
in T
Fig.
Clea



Superior Performance and Proven Biocompatibility


Comparative Evaluation of Clearance

Clearances of new Toraysulfone TS series and BS series (current product) were compared in 8 end stage renal disease (ESRD) patients (crossover study).


Clearances of phosphate and α 1-MG in TS series were significantly higher than BS series. (Fig. 6)

Comparison of clearances in 6 ESRD patients (crossover) between TS series and another polysulfone dialyzer (APS), clearances in uric acid and α 1-MG were also significantly higher with Toraysulfone TS series. (Fig. 7)

Comparison of Clearance in TS and BS



Comparison of Clearance in TS and APS



Comparison of Biocompatibility

Significant change in C3a and leukocyte counts were observed with FB (CTA membrane). In contrast, less significant changes were observed with Toraysulfone membrane. This suggests that the Toraysulfone membrane has better biocompatibility than that of FB. (Fig. 8, 9)

Fig. 8: Changes in C3a comparing Toraysulfone and FB





Fig. 9: Changes in leukocyte counts comparing Toraysulfone and FB




Technical Specifications and Data

Type	TS-1.3S	TS-1.6SL	TS-1.8SL	TS-2.1SL	TS-1.3U	TS-1.6UL	TS-1.8UL	TS-2.1UL
Housing Material	Polycarbonate							
Fibers Material	Polysulfone							
Inner diameter (μm)	200							
Membrane thickness (μm)	40							
Effective surface area (m ²)	1.3	1.6	1.8	2.1	1.3	1.6	1.8	2.1
Potting Material	Polyurethane							
Sterilization	Gamma-ray Irradiation							
Blood Volume (mL)	84	95	105	124	85	95	108	125
Clearance in vitro (mL/min)*								
Urea	193	195	198	198	193	196	198	199
Creatinine	187	193	195	195	185	192	196	197
Phosphate	179	192	194	195	180	193	196	196
Vitamin B ₁₂	140	156	164	168	140	162	167	171
Inulin	104	124	129	138	110	131	140	142
UFR in vitro (mL/hr, at 13.3kPa (100mmHg))**	4,400	4,900	5,000	5,200	4,300	4,900	5,100	5,500
Max. TMP (kPa (mmHg))	66 (500)							

* Clearances are measured with aqueous solution. Q_B: 200 ± 4mL/min, Q_D: 500 ± 10mL/min, Q_F: 10mL/min, Temp.: 37 ± 1°C
 ** UFRs are measured data with bovine blood. (Ht 30 ± 2%, TP 6g/dL) Q_B: 200mL/min, TMP: 13.3kPa (100mmHg), Temp.: 37°C
 Allowable ranges
 Blood volume; ±13%. Clearances, Urea & Creatinine; ±6%, Phosphate & Vitamin B₁₂; ±13%, UFR in vitro; ±15%
 "Instructions for Use" should be read thoroughly prior to the use of these medical devices.
 Specifications and designs are subject to change without notice for improvements.


Clearance



TS-1.3S	TS-1.3U		
Q _B	200 300 400 500	Q _B	200 300 400 500
UR; Urea	193 255 289 320	UR; Urea	193 257 294 323
CR; Creatinine	187 235 260 285	CR; Creatinine	185 238 264 285
PO ₄ ; Phosphate	179 223 250 268	PO ₄ ; Phosphate	180 224 250 268
VB ₁₂ ; Vitamin B ₁₂	140 161 172 187	VB ₁₂ ; Vitamin B ₁₂	140 164 178 187
IN; Inulin	104 117 131 140	IN; Inulin	110 123 135 143
TS-1.6SL	TS-1.6UL		
Q _B	200 300 400 500	Q _B	200 300 400 500
UR; Urea	195 271 320 348	UR; Urea	196 274 326 358
CR; Creatinine	193 255 293 315	CR; Creatinine	192 259 301 327
PO ₄ ; Phosphate	192 248 283 307	PO ₄ ; Phosphate	193 250 287 311
VB ₁₂ ; Vitamin B ₁₂	156 187 207 221	VB ₁₂ ; Vitamin B ₁₂	162 194 216 232
IN; Inulin	124 142 156 175	IN; Inulin	131 152 170 184
TS-1.8SL	TS-1.8UL		
Q _B	200 300 400 500	Q _B	200 300 400 500
UR; Urea	198 275 332 363	UR; Urea	198 277 332 367
CR; Creatinine	195 262 306 333	CR; Creatinine	196 264 308 340
PO ₄ ; Phosphate	194 254 292 317	PO ₄ ; Phosphate	196 258 297 324
VB ₁₂ ; Vitamin B ₁₂	164 197 218 234	VB ₁₂ ; Vitamin B ₁₂	167 202 226 248
IN; Inulin	129 153 165 182	IN; Inulin	140 162 182 202
TS-2.1SL	TS-2.1UL		
Q _B	200 300 400 500	Q _B	200 300 400 500
UR; Urea	198 279 335 369	UR; Urea	199 280 336 371
CR; Creatinine	195 268 313 342	CR; Creatinine	197 268 314 344
PO ₄ ; Phosphate	195 258 296 321	PO ₄ ; Phosphate	196 259 299 326
VB ₁₂ ; Vitamin B ₁₂	168 204 230 249	VB ₁₂ ; Vitamin B ₁₂	171 207 235 253
IN; Inulin	138 160 180 196	IN; Inulin	142 166 185 204

Measured data with aqueous solution

UFR



References

- 1) Sakai, K., taken from the presentation made at the 48th Annual Meeting of the Japanese Society for Dialysis Therapy, Jun. 2003 (Fig.1, Table 1)
- 2) Hayama, M., Yamamoto, K., Kohori, F., Sakai K., et al., Biomaterials 25 (2004) 1019-1028 (Fig.1, Table 1)
- 3) Sakai, Y., et al., taken from the presentation made at the 22nd Annual Meeting of ISBP, Sep. 2004 (Fig. 5)
- 4) Miyaji, H., et al., taken from the presentation made at the 20th Annual Meeting of HPM, Mar. 2004 (Fig. 6, 7)
- 5) Onishi, N., et al., taken from the presentation made at the 31st Shikoku Dialysis Therapy Research Group, 1997 (Fig. 8, 9)